Lynx™ Repeater

ENGLISH: Repeater and receiver of doorbell for the Lynx system.....	page 2
DANSK: Forstærker og modtager af dørsignal til Lynx-systemet.....	side 6
NEDERLANDS: Repeater en ontvanger voor deurbel van Lynx-systeem.....	pag. 10
NORSK: Repeater og mottaker av dørsignal til Lynx-systemet.....	side 14
SVENSKA: Repeater och mottagare av dörrsignal till Lynx-systemet	sid 18


Ordering number: 200 21 06


Repeater for the Lynx system

The repeater forms part of Lynx – our wireless alert system for people with a hearing impairment. Lynx uses various types of transmitters and receivers that alert you to noises such as a ringing telephone or doorbell, or a fire alarm.


The repeater has *two* functions. *Firstly*, it functions as a **repeater**, i.e. it receives, amplifies and transmits the signal onwards so that it reaches further. With one or several Repeaters you enlarge the range and security of the system.

Secondly, the Repeater can function as an **acoustic door alert**, which uses a powerful signal to alert the user.

Overview

Installing the Repeater

1. *Choosing location:*

Select a suitable position for the Repeater in a central location in the home.

It should always be installed on a flat surface in the vicinity of an electric socket, ideally on the wall approx. 150-180 cm above the socket. Avoid positioning in corners and in the vicinity of metal cabinets and similar materials as it can prevent the wireless transmission of the signals.

2. *Channel setting:*

Check that the Repeater and all transmitters and receivers in the system have the same channel setting (see next page).

3. *Connection to the mains:*

Connect the mains adapter cable to the socket under the Repeater. Put the cord in the electric socket.

4. *Function test:*

Install the transmitters that you are going to use. **Read the user manual for the transmitters that you are installing.**

Test the system by activating the transmitters in the system one at a time and check that the receiver is displaying the correct alarm and that it can be reached in all places.

5. *Volume:*

Set the required volume for the door alert (strong, weak or off) with the volume control.

6. *Ready for use:*

The Repeater does not have an on/off switch, but is now ready for use.

Changing channels

The transmitters in the Lynx system transmit the signal wirelessly to the Repeater, which transmits the signal onward. They can transmit on a total of 16 different channels (0-9, A-F). The Repeater and all transmitters and receivers must be set to the same channel so that they can communicate. They are preset to channel 9, and then operate together without alteration.

In some cases, it may be necessary to change the channel setting. If two or more users live near one another, the alarms can interfere. This can be rectified by changing the channel settings so that they are different for each person. Any channel can be selected except channel 0. (Channel 0 is reserved for a common alert such as a fire alarm.)

If you need to change the channel on the Repeater, proceed as follows:

The channel selector is on the rear. Carefully screw on the channel selector using a screwdriver.

To change channels on the receiver, see user manual.

Troubleshooting

Start by reading the user manual. Check that the connection is correct and that all transmitters and receivers in the system are set to the same channel. Then look in the table below.

Type of fault	Possible cause	Action
The receiver does not receive an alert from one of the transmitters, but the others are functioning properly.	<ol style="list-style-type: none"> 1. Too great a distance between the transmitter and the Repeater. 2. The transmitter has a flat battery. 3. The transmitter is defective. 	<ol style="list-style-type: none"> 1. Move them closer and check whether the signal is received. 2. Change the battery (see transmitter user manual). 3. Contact the place of purchase for repair.
The receiver does not always receive alerts from some transmitters.	There is too great a distance between the Repeater and the receiver or between the transmitters and the Repeater.	Try moving the Repeater to another place or increase the number of Repeaters.
The receiver often receives false alarms.	Someone in the vicinity has the same system installed.	Change to another channel (see above).
No sound from the speaker.	<ol style="list-style-type: none"> 1. The volume control is set at zero. 2. The Repeater only gives an acoustic alarm for door alert (other alarms give no sound but the Repeater is still functioning). 	<ol style="list-style-type: none"> 1. Increase the volume. 2. None.

Cleaning and maintenance

Vacuum the outside of the Repeater with a soft brush. Use a damp cloth to wipe it on the outside. Water must not enter inside the cover. Never use cleaning products or solvents.

Repairs

If your product malfunctions, it must be repaired by a qualified technician. Do not attempt to open the case of the device since this would invalidate the warranty. If your product requires service, please contact your hearing care professional for assistance.

Waste disposal

This product contains electrical or electronic equipment and should be disposed of carefully in the interests of your safety and the environment. Please contact your local hearing care professional concerning disposal of the product.

General precautions

- Different types of radiation, e.g. from NMR or CT scanners, may damage the device. Therefore, do not wear the device in the vicinity of this type of equipment. Other types of radiation (burglary alarms, room surveillance systems, radio equipment, mobile telephones, etc) will not damage the device. These types of radiation could, however, momentarily affect the sound quality or create noise in the device.
- **WARNING:** External devices connected to the product must be safe according the requirements of EN 60601-1-1 or EN 60065, as appropriate.

Technical data

Frequency:	433.92 MHz.
Modulation:	AM pulsing.
Battery:	1 x 9 V type 6LR61 or IEC 6F22. Provides approx. 16 hours back-up if the Repeater has been connected for at least one week.
Power consumption:	<i>Standby:</i> approx. 0.2 mA. <i>Transmission:</i> approx. 16 mA.
Range:	approx. 100 m in an open field.
Channels:	16.
Dimensions:	170 x 70 x 45 mm.
Weight:	approx. 200 g.

Any issues relating to the EU Medical Device Directive 93/42/EEC should be directed to GN ReSound A/S.

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INTERNATIONAL INQUIRIES:**GN ReSound AB ALD Division**
6930 AB, Westervoort
SVERIGE:*Postadress:***GN ReSound AB**
*Service:***GN ReSound AB ALD Division**
MANUFACTURER:*in the sense of EU Medical**Device Directive 93/42/EEC:***GN ReSound A/S**