Lynx™ TX-1 Transmitter


Ordering number: 200 30 24

ENGLISH: Universal transmitter to our alerting system Lynx.....	page 3
DEUTSCH: Universal-Sender für das Lynx Alarmsystem	Seite 8
FRANÇAIS: Emetteur universel vers le système d’alerte Lynx	page 13
ITALIANO: Trasmettitore universale per il nostro sistema di allarme Lynx	Pag. 18
NEDERLANDS: Universele zender voor het Lynx-systeem.....	Pag. 23
NORSK: Universalsender til varslingsystemet Lynx.....	side 28
SVENSKA: Universalsändare till varseblivningssystemet Lynx	sid 33

Universal transmitter for the Lynx system

Lynx TX-1 forms part of Lynx – our wireless alert system for people with impaired hearing. The Lynx uses various transmitters and receivers that alert you to noises such as a ringing telephone or doorbell, or the first signs of a fire.

Lynx TX-1 is a universal transmitter that can detect and transmit many different types of signals such as a telephone signal, doorbell, fire alarm, activated photocells and activated magnetic switches. Signals are sent wirelessly to the receiver. One or more repeaters may also be installed in order to amplify the alarm signal and thus extend its range.


Connection

There are two different connection options for Lynx TX-1. One type of connection is for the telephone (via a modular switch) and the other is the connection of a signal cable to the signal source (via a connection block).

Lynx TX-1 and telephone

The accessories supplied are a cable with modular switches and an intermediate plug.

1. Connect the cable to the modular switches between Lynx TX-1 and the intermediate plug.
2. Insert the intermediate plug in the telephone socket.
3. Then connect the telephone plug to the intermediate plug.


Lynx TX-1 and other equipment

Other types of signal sources can be connected via the three connection blocks under the cover. Use the cable supplied for connection.

If your signal source is a short-circuiting signal (eg doorbell) it is connected via connection blocks 2 and 3.

Equipment which supplies a voltage when it is activated can also be connected. It is connected between connection blocks 1 and 2. The voltage activating Lynx TX-1 must be between 0.7-24 V, either alternating current or DC voltage.

Fitting Lynx TX-1

Screws are supplied in the package which are used to fit Lynx TX-1 onto the wall. Lynx TX-1 can also be attached using the Velcro tape supplied. If Lynx TX-1 is used at different locations in the house and you want to be able to move it rapidly, Velcro tape can be attached at different locations.

When the sound source is activated

When the sound source is activated, a wireless signal is sent to the receiver installed in the Lynx system. The alert is given by vibration and symbols, loud noise or flashing light, depending on the type of receiver.

Read the user manual for the receiver installed!

Flat battery

Lynx TX-1 is powered by a battery with a normal service life of approximately 5-10 years. When battery power begins to fail, Lynx TX-1 will send out an alert. It sends a flat-battery warning along with the alarm to the receiver in the alerting system.

Read the user manual for the receiver installed!

Changing the battery

1. *Open the blue cover.*

2. *Change the battery.*

Setting the alarm symbol and changing the channel setting

Standard settings:

Test switch for installation. When it is set to ON, Lynx TX-1 transmits continuously. This makes it easier to find the most suitable location for the receiver.

Channel selector switch Alarm selector switch

Lynx TX-1 is equipped with a transmitter that sends signals wirelessly to the receiver installed in the Lynx system. The transmitter is located under the blue cover and has an 8-pole microswitch that is used for two purposes:

- **Alarm symbol** (fire, telephone, door, etc). Lynx TX-1 is pre-set for door but another suitable symbol can be selected if required. The symbols that can be selected depend on the receiver that is being used. If you wish to change the type of alarm, set switches 6, 7 and 8 on Lynx TX-1 as shown in the table alongside.
- **Channels.** There is a total of 16 different channels (0-9, A-F). Transmitters and receivers in the system must be set to the same channel so that they can communicate. Channel 0 is reserved for a common fire alarm and should not be used. Lynx TX-1 and all receivers in the Lynx system are all pre-set to channel 9, and then act together without alteration.

Changing the channel setting

In some cases, it may be necessary to change the channel setting. If two or more users live near one another, the alarms can interfere. This can be rectified by changing the channel settings so that they are different for everyone. Any channel can be selected except Channel 0. If you wish to change channel, set the switches on Lynx TX-1 as shown in the table alongside. To change channels on the receiver, see user manual.

Setting the alarm selector microswitches

	Switch:		
	6	7	8
Telephone	On	On	Off
Door	On	Off	On
Fire	Off	On	On
Baby	Off	Off	Off
Pager	On	Off	Off
Optional	Off	Off	On
Public alert	Off	On	Off

Setting the channel selector microswitches

Channel:	Switch:				
	1	2	3	4	5
0	Off	Off	Off	Off	Off
1	On	Off	Off	Off	Off
2	Off	On	Off	Off	Off
3	On	On	Off	Off	Off
4	Off	Off	Off	On	Off
5	On	Off	Off	On	Off
6	Off	On	Off	On	Off
7	On	On	Off	On	Off
8	Off	Off	Off	Off	On
9	On	Off	Off	Off	On
A	Off	On	Off	Off	On
B	On	On	Off	Off	On
C	Off	Off	Off	On	On
D	On	Off	Off	On	On
E	Off	On	Off	On	On
F	On	On	Off	On	On

Common channel for fire alarm

Pre-set channel on Lynx TX-1 and all receivers

Troubleshooting

Type of fault	Action:
Lynx TX-1 fails to send out alarm.	<ol style="list-style-type: none"> 1. Check that the battery is OK. 2. Check that the cables are properly fitted.
Receiver fails to respond.	Check the range. Try to reduce the distance between Lynx TX-1 and the receiver or install a Repeater to extend the range.

Cleaning and maintenance

Vacuum outside of Lynx TX-1 with a soft brush several times per year. Use a lightly moistened cloth to dry off the exterior. Never use cleaning products or solvents.

Repairs

If the product malfunctions it must be repaired by an authorised technician. Do not attempt to open the case of the device since this would invalidate the warranty. If your product requires service, please contact your hearing care professional for assistance.

Waste disposal

This product contains electrical or electronic equipment and should be disposed of carefully in the interest of your safety and the environment. Please contact your local hearing care professional concerning disposal of the product.

General Precautions

- Low-energy radiation (eg from burglary alarms, room surveillance systems, radio equipment, mobile telephones) will not damage the device.
- Warning: External devices connected to the product must be safe according the requirements of EN 60601-1-1, EN 60335-1 or EN 60065, as appropriate.

Technical data

Modulation:	AM.
Frequency:	433.92 MHz.
Range:	>50 m free line of sight to repeater or receiver. Shorter range indoors.
Channels:	16 via microswitches.
Functions:	1 of 7 options is available.
Telephone input:	Modular jack.
Other inputs:	0.7-24 V AC/DC and short circuit.
Battery:	1 x 9 V type 6LR61 or IEC 6F22.

Any issues relating to the EU Medical Device Directive 93/42/EEC should be directed to GN ReSound A/S.

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	Schakelaar:		
	6	7	8
Telefoon	On	On	Off
Deur	On	Off	On
Brand	Off	On	On
Baby	Off	Off	Off
Pieper	On	Off	Off
Optioneel	Off	Off	On
Openbaar alarm	Off	On	Off

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Kanaal:	Schakelaar:				
	1	2	3	4	5
0	Off	Off	Off	Off	Off
1	On	Off	Off	Off	Off
2	Off	On	Off	Off	Off
3	On	On	Off	Off	Off
4	Off	Off	Off	On	Off
5	On	Off	Off	On	Off
6	Off	On	Off	On	Off
7	On	On	Off	On	Off
8	Off	Off	Off	Off	On
9	On	Off	Off	Off	On
A	Off	On	Off	Off	On
B	On	On	Off	Off	On
C	Off	Off	Off	On	On
D	On	Off	Off	On	On
E	Off	On	Off	On	On
F	On	On	Off	On	On

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	Bryter:		
	6	7	8
Telefon	On	On	Off
Dør	On	Off	On
Brann	Off	On	On
Baby	Off	Off	Off
Pager	On	Off	Off
Valgfri	Off	Off	On
Publikumsalarm	Off	On	Off

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Kanal:	Bryter:				
	1	2	3	4	5
0	Off	Off	Off	Off	Off
1	On	Off	Off	Off	Off
2	Off	On	Off	Off	Off
3	On	On	Off	Off	Off
4	Off	Off	Off	On	Off
5	On	Off	Off	On	Off
6	Off	On	Off	On	Off
7	On	On	Off	On	Off
8	Off	Off	Off	Off	On
9	On	Off	Off	Off	On
A	Off	On	Off	Off	On
B	On	On	Off	Off	On
C	Off	Off	Off	On	On
D	On	Off	Off	On	On
E	Off	On	Off	On	On
F	On	On	Off	On	On

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	Brytare:		
	6	7	8
Telefon	On	On	Off
Dörr	On	Off	On
Brand	Off	On	On
Baby	Off	Off	Off
Anhörlarm	On	Off	Off
Valfri	Off	Off	On
Publikt larm	Off	On	Off

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Kanal:	Brytare:				
	1	2	3	4	5
0	Off	Off	Off	Off	Off
1	On	Off	Off	Off	Off
2	Off	On	Off	Off	Off
3	On	On	Off	Off	Off
4	Off	Off	Off	On	Off
5	On	Off	Off	On	Off
6	Off	On	Off	On	Off
7	On	On	Off	On	Off
8	Off	Off	Off	Off	On
9	On	Off	Off	Off	On
A	Off	On	Off	Off	On
B	On	On	Off	Off	On
C	Off	Off	Off	On	On
D	On	Off	Off	On	On
E	Off	On	Off	On	On
F	On	On	Off	On	On

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INTERNATIONAL INQUIRIES:**GN Resound AB ALD Division**
GN Resound A/S