Transett 32 (all versions)

USER MANUAL FOR:	ART. NO:	ENGLISH	page 2
Transett 32	200 10 18	SVENSKA	sidan 16
Transett 32 24V	200 10 22	NORSK	side 30
Transett 32 Phase Shift 90°	200 10 21	NEDERLANDS	Pag. 44
Transett 32 Phase Shift 180°	200 10 26		
Transett 32 Phase Shift Card	877 46 31 00		

Contents

1. Three versions for different purposes.....	3
2. Mounting the loop amplifier	3
Connection terminal.....	4
3. Loop installation	4
Choosing loop wire	4
Single loop installation.....	5
Divided loop installation.....	6
Multiple loops installation	7
(only Transett 32 Phase Shift 90° & 180°)	
4. Connecting and adjusting sound signals	9
AGC	9
Adjusting loop current.....	11
5. Installing a Phase Shift Card in Transett 32	12
6. Technical data.....	14
7. Troubleshooting	15
8. Maintenance and safety	15

1. Three versions for different purposes

Transsett 32 is a loop amplifier with robust design specifically suitable for use in churches, meeting halls and schools. There are three versions available:

- **Transsett 32** is our standard version and used for single loop implementations. Transsett 32 can be upgraded to a Transsett 32 Phase Shift 90° or 180° using a **Transsett 32 Phase Shift Card**.
- **Transsett 32 Phase Shift 90°** and **Transsett 32 Phase Shift 180°** are two loop amplifiers that combined generate two fields with 90 degree phase shift between them. This creates a homogenous field in the room with strong field attenuation outside the loop area.
- **Transsett 32 24V** is specifically designed for 24 V environments in transport vehicles (e.g. train, bus and boat).


This user manual covers all versions. If not otherwise specified, the information is applicable to all versions.

2. Mounting the loop amplifier

The loop amplifier can be placed on a shelf or mounted on the wall. It is equipped with rubber feet to ensure safe placement on a shelf. The package includes screws and plugs for wall mounting in an arbitrary orientation. It is normal for the loop amplifier to run slightly warm. It is therefore recommended to place it in an open location and it must not be covered by any object to ensure best performance. The amplifier contains a thermal device which automatically reduces the signal level if over-heated, thus supplying a weaker signal to the hearing instrument.

Connection terminal

The loop amplifier can be rapidly and easily fitted and removed, for example for servicing. The signal sources are connected to a detachable connection strip, which can easily be removed and attached. (For connection of signal sources, see page 10, step 5).

3. Loop installation

Choosing loop wire

All Transsett 32 versions can be used for loops up to approximately 180 m². For optimal result the loop, wire and amplifier parameters should be well matched. The resistance in the loop should be approximately 1.5 Ω regardless of the loop area. Recommended field strength is 100 mA/m 1.2 m above the floor (which is normal listening height). Use the information below for a rough estimation. To obtain exact values, we recommend using a field strength meter.

Area covered (m ²)	Wire dimension (mm ²)	Transsett 32 potentiometer position
180	1.5	180
90	0.4-0.75	90
45	0.4	45

Set the potentiometer value according to the table above, using a screw driver.


Single loop installation

Using a single loop around the room, placed close to the floor, is the most common type of loop installations.

In most cases a single loop installation provides a satisfactory result, but it does have a few disadvantages:

- A low attenuation of the field outside of the loop.
- Above the loop wire the field is only horizontal.
- When large areas are covered, the field may be weaker in the center.

Divided loop installation


In environments where sharper field attenuation outside of the loop is desired, a divided loop is a better alternative.

Advantages:

- Creates a more homogenous field within respective loop area.
- The field is better attenuated outside of the loop if the parts are of equal size and the fields have reversed directions.

Disadvantage:

- At the borderline between the two loops a cancellation effect occurs between the two loops, i.e. no signal is received.

Multiple loops installation (only Transett 32 Phase Shift 90° & 180°)

When sharp field attenuation outside the loop and a regular field within is desired, the best solution is multiple loops fed with phase shifted currents. Phase shifted currents can be obtained by feeding each of the two wire loops by two separate loop amplifiers: Transett 32 Phase Shift 90° & 180°. Each of the two wire loops may in turn be configured as multiple sub-loops of equal size (always a multiple of two i.e. 2, 4, 8, 16 loops etc). These sub-loops should be placed edge-to-edge. The more loops used, the better the attenuation of the field outside the loops. The different loop areas shall be of the same size for maximum field attenuation outside the loop area.

Advantages:


- Strong field attenuation outside the loop area.
- Very homogenous field within the loop.
- The hearing instrument works well in the whole area covered, independent of its orientation.

Connecting Transett 32 Phase Shift 90° & 180°

- It does not matter which one of Transett 32 Phase Shift 90° or 180° you use for signal input, but the one used for signal input should via Line Out feed the other one.
- Anyone of the 40 mV inputs can be used (i.e. Line 1, Line 2 or Line 3).
- One of the input terminals must be grounded to 8.

Example

An installation using four loops. In this installation connection to "Line 3" is used. The numbers in this figure refer to the numbers on the connection terminal.

4. Connecting and adjusting sound signals

AGC

Transett 32 (all versions) are equipped with Automatic Gain Control (AGC). This is to ensure that the signal input from each source is within an acceptable level and to prevent the signal to the loop from being too strong. AGC characteristics can be set to three values: Slow, Fast & Off. The default setting is *Fast*, which is recommended for most listening situations.

Single loop installation using one Transett 32 or one Transett 32 24V:
Set desired AGC and adjust loop current according to the description on pages 9-12.

Multiple loop installation using both Transett 32 Phase Shift 90° & 180°:
First, set desired AGC and adjust loop current on Transett 32 Phase Shift 90° according to the description on pages 9-12. Then, repeat the procedure for Transett 32 Phase Shift 180°.


Setting AGC:

Recommended value: Fast (default).


1. Make sure that the loop amplifier is NOT connected to a mains power outlet. Then, unscrew the hood and remove it.
2. Locate the AGC jumper on the IC-holder.

Transett 32 & Transett 32 24V

Transett 32 Phase Shift 90° & 180°

3. Set the AGC characteristics to the appropriate value (see table below).


AGC Characteristic	Jumper position
Slow (e.g. music)	9 & 10
Fast (e.g. speech, music) default setting	11 & 12
Off (when using an even signal source, e.g. a CD)	10 & 11

4. Replace the hood.

5. Connect the desired signal source to the input with the most appropriate sensitivity.

All inputs are electronically balanced. When using a *balanced signal source*, connect it with a shielded pair cable to the connection terminal of the loop amplifier (numbers 1 and 2 when using Line 2). Connect the earth cable to the ground terminal (number 8 on the connection terminal).

When using an *unbalanced signal source*, connect it with a shielded single threaded cable to the connection terminal of the loop amplifier (numbers 1 or 2 when using Line 2). Connect the earth cable to the ground terminal (number 8 on the connection terminal).

6. Connect the loop amplifier to a mains power outlet.

There is no mains power switch on the loop amplifier. It can be left on round-the-clock. Power consumption is very limited, especially when idle, i.e. in the absence of a programme signal.

7. Add sound. The sound source should be the one you are planning on using in the real situation. If it is a TV, then use the TV as source. If it is a microphone, then use the same microphone and talk into it etc.

8. Turn the volume control for the appropriate input slowly clockwise until the yellow AGC LED lights up. Then turn the control slightly counter-clockwise until the yellow light goes out. AGC should flash at times but never be lit continuously.

9. Repeat the procedure for all signal sources.

Adjusting loop current

1. For a rough estimation of the potentiometer position, see table on page 4. Then turn the volume control "LOOP AREA" clockwise until the green lamp "LOOP CURRENT" is lit up indicating that a suitable current flows through the loop.

2. If the red lamp “OVERLOAD” starts to flash, the amplifier is overloaded. If so, turn back volume control “LOOP AREA” slightly.
3. For fine adjustment we recommend using a field strength meter to measure the loop field intensity. We recommend an average field intensity value (measured at 1.2 m above the floor) of 100 mA/m and maximum 400 mA/m in peaks. Adjust, if necessary, with the volume control “LOOP AREA” as described above.

When fine adjusting the field strength in installations with both Transsett 32 Phase Shift 90° & 180°, make sure that you measure the field strength in the middle of the loop connected to the amplifier you are adjusting.


5. Installing a Phase Shift Card in Transsett 32

Transsett 32 can be upgraded to a Transsett 32 Phase Shift 90° or 180° using a Transsett 32 Phase Shift Card.


1. Make sure that Transsett 32 is NOT connected to a mains power outlet. Then, unscrew the hood and remove it.


2. Locate the IC-holder and remove the two jumpers from it.


3. Gently press the Phase Shift Card into the same IC-holder.

4. Fasten the Phase Shift Card with a screw.

5. Connect the jumper to obtain a 90 or 180 degree phase shift.
Pin 1-1: 90°
Pin 2-2: 180°


6. Replace the hood.

6. Technical data

Power supply (32 and Phase Shift):	210-264 V 50-60 Hz, nominally 230 V.
Power supply (32 24V):	20-30 V DC.
Inputs (DIN-connector):	Microphone 0.5 mV-15 mV/8 kΩ. Balanced with phantom 15 V. Line 1 level 40 mV-1V/20 kΩ. Electronically balanced.
Inputs (connection terminal):	Line 2, and 3 level 40 mV-1V/20 kΩ. Electronically balanced. Line 4: level 400 mV-10 V/ 45 kΩ. Electronically balanced.
Outputs:	Loop output 30 VA with 1.5 Ω loop resistance. 24 V DC, 300 mA for accessories. Line out 600 mV. Power supply for door chime and buzz indicator on loop.
Distortion:	1% at 30 VA.
Frequency range:	100-9000 Hz, - 3dB rel. 1000 Hz. Roll off > 6dB/octave above 9000 Hz.
Signal/noise ratio:	>50 dB for microphone input. >60 dB for line inputs.
Treble boost:	+5 dB at 4000 Hz rel. 1000 Hz cut-off point at 2000 Hz.
Bass boost:	+5 dB at 200 Hz rel. 1000 Hz cut-off point at 500 Hz.
AGC:	40 dB dynamic range (used to maintain a constant average current in the loop). Attack-time 30 ms. Decay-time 30 seconds, adjustable with jumper.
Indicators:	Yellow LED lit on activation of AGC. Green LED indicates current loop. Red LED indicates imminent overload.
Overload protection:	Thermal; reduce drive level.
Power consumption:	70 W max. Idle state: <5 W.
Dimensions:	258 x 120 x 75 mm.
Weight:	2.5 kg / 0.8 kg (24V).
Colour:	Grey.
Ambient temperature:	10-40 °C.
All versions of Transsett 32 meet all the requirements of IEC60118-4.	

CE Any issues relating to the EU Medical Device Directive 93/42/EEC should be directed to GN ReSound A/S.

7. Troubleshooting

Type of fault	Cause	Action
The sound is faint.	The loop field is too weak.	<ol style="list-style-type: none"> 1. First check that the volume control has been set correctly as described in steps 8-9 on page 11. 2. Then check that the loop field has been adjusted as described in steps 1-3 on page 11-12.
The sound is barely audible.	The loops “cut out” one another (for phase shift/un-correlated loops).	The phase offset difference must be 90 degrees (see page 7 and also page 13 step 5 to check the jumper settings).
The sound is too loud or distorted.	The loop field is too strong.	<ol style="list-style-type: none"> 1. First check that the volume control has been set correctly as described in steps 8-9 on page 11. 2. Then check that the loop field has been adjusted as described in steps 1-3 on page 11-12.
The sound fluctuates unnaturally with music.	AGC is set to the “fast” position.	Set AGC to the “slow” or “off” position (see page 10, step 3.)
The sound fades out but comes back slowly during speech.	AGC is set to the “slow” position.	Set AGC to the “fast” or “off” position (see page 10, step 3.)
Transset 32 is “completely dead”.	One or both fuses may have blown.	Check the fuses: 1) Make sure that Transset 32 is not connected to 220 volts! 2) Unscrew the casing as described in step 1, page 9. 3) The fuses are under the card (on the same side as the IC holder for AGC, see step 2, page 9).

8. Maintenance and safety

Cleaning and maintenance

Use a vacuum cleaner with a soft brush to clean the product on the outside. Use a damp cloth to wipe it on the outside. Never use cleaning products or solvents.

From time to time check that the cables and connections are intact.

Repairs

If your product malfunctions, it must be repaired by a qualified technician. Do not attempt to open the case of the device since this would invalidate the warranty. If your product requires service, please contact your hearing care professional for assistance.

Waste disposal

This product contains electrical or electronic equipment and should be disposed of carefully in the interest of your safety and the environment. Please contact your local hearing care professional concerning disposal of the product.

General precautions

- Low-energy radiation (eg from burglar alarms, room surveillance systems, radio equipment, mobile telephones) will not damage the device. They could, however, momentarily affect the sound quality or create noise in the device.
- **Warning:** External devices connected to the product must be safe according the requirements of EN 60601-1-1, EN 60335-1 or EN 60065, as appropriate.

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SVERIGE:*Postadress:***GN ReSound AB**
*Service:***GN ReSound AB ALD Division**
INTERNATIONAL INQUIRIES:**GN ReSound AB ALD Division**
MANUFACTURER:*in the sense of EU Medical Device**Directive 93/42/EEC:***GN ReSound A/S**